

Protocol: Miniprep Plasmid DNA Extraction

This is a generic protocol for plasmid DNA extraction and purification. Many plasmid miniprep kits are sold that operate on similar principles. You should follow the specific instructions for the kit that you buy.

Material	Amount	Material	Amount
Resuspension Buffer	250 μ l	Miniprep Spin Column	1
Cell Lysis Buffer	250 μ l	Ethanol Wash Solution	1 mL
Neutralization Buffer	550 μ l	Elution Buffer or Water	50 μ l

The night before:

- 1) Start a 5 mL overnight culture of the bacteria carrying your plasmid.**
Include a selective antibiotic.
- 2) Centrifuge 2 mL of culture for 2 minutes at 6000 rpm.**
Keep the cell pellet and discard the supernatant.
- 3) Resuspend the cell pellet in 200 μ l Resuspension Buffer.**
- 4) Add 250 μ l Cell Lysis Buffer.**
Mix gently and wait 5 minutes for lysis to proceed.
- 5) Add 350 μ l Neutralization Buffer.**
Mix gently.
- 6) Centrifuge the lysed cells for 5 minutes at maximum speed.**
Keep the supernatant and discard the cell debris.
- 7) Add the supernatant to a Miniprep Spin Column and spin for 1 minute.**
Keep the column and discard the flow-through.
- 8) Add 500 μ l Ethanol Wash Solution to the column and spin for 1 minute.**
Keep the column and discard the flow-through.
- 8) Again add 500 μ l Ethanol Wash Solution to the column and spin for 1 minute.**
Keep the column and discard the flow-through.
- 9) Spin the column for 2 minutes at maximum speed to dry it.**
Ensure that no ethanol remains on the column
- 10) Transfer the column to a clean Eppendorf tube.**
- 11) Add 30-100 μ l of water directly to the column membrane.**
Wait a few minutes to allow the DNA to elute.
- 12) Spin the column at max speed for 1 minute to collect the plasmid DNA.**
Keep the flow-through and discard the column.